

Transient Characteristics of GDI Injectors using Statistical Extinction Tomography




Yudaya Sivathanu, Jongmook Lim, and Varun Kulkarni

En'Urga Inc.

1201 Cumberland Ave., West Lafayette, IN 47906


This work was funded by the National Science Foundation and Delphi Corp., USA. The assistance of Mr. Vinay Prasad (Woodward Inc.) with the experiments is also acknowledged.





Outline


- **Background**
- **Experimental methods**
- **Results**



Background


- *Spray characteristics affect engine performance (Speigel and Spicher, 1992, Davy et al., 2000, Zhoa et al., 2002)*
- *Multi-hole configuration better than pressure swirl atomizers (Fujieda et al., 1995; Arcoumanis et al., 1998)*
- *Spray characteristics vary with ambient pressure and fuel temperature (Schmitz et al., 2002; Parish et al., 2010)*

Objective of this study is to investigate the use of extinction tomography to characterize dense multiple plume sprays



Experimental Methods

Optical Configuration




- Single axis extinction measurement at 9.4 KHz
- Planar extinction through 512 parallel paths
- Injector rotated 12 times to obtain multiple view data

Experimental Arrangement



Setscan AP400 patternator


- *250 mm nominal ID*
- *100 mm silica windows*
- *~ 1m high vessel to reduce spray bounce*
- *Injector rotated by 15 deg for a total of 12 views*
- *Inert nitrogen atmosphere*
- *Fuel heated upto 90 °C*
- *5 injections per angle*



Principle of Operation

- Path integrated extinction of laser sheets
- 12 view angles, 15 degrees apart
- 5 injection events ensemble averaged at each view angle
- 512 parallel path measurements at each view angle
- Local extinction coefficients in the plane obtained by deconvolution (MLE – Vardi and Lee, 1993)
- For dielectric liquids, local extinction coefficient is proportional to the **drop surface areas per unit volume**

Importance of Surface Area Density



Control volume

Drop size varied from 1 to 100 micron
Velocity varied from 1 to 10 m/s
Number of drops varied from 1E5 to 1E7
Calculate evaporation rate

For static drop $\dot{m} = \frac{4\pi k_g r_s}{C_{pg}} \ln(B+1)$

For dynamic drop $\dot{m}_{\text{dynamic}} = \frac{N_u}{2} \dot{m}_{\text{static}}$


$$B = \frac{C_{pg}(T_\alpha - T_b)}{h_{fg}}$$

$$C_{pg} = C_{pF}(\bar{T}) \quad \text{and} \quad k_g = 0.4k_F(\bar{T}) + 0.6k_\infty(\bar{T}); \quad \bar{T} = (T_b + T_\infty)/2$$

$$N_u = 2 + \frac{0.555 \text{Re}^{1/2} \text{Pr}^{1/3}}{\left[1 + 1.232/(\text{RePr})^{4/3}\right]^{1/2}}$$

*An Introduction to Combustion: Concepts and Applications,” Stephen R. Turns, Mc. Graw Hill, 2000

Importance of Surface Area Density



Correlation coefficient (R) of spray parameters with evaporation

Mass flux $R = 0.903$

Velocity $R = -0.239$

Diameter $R = 0.681$

Surface area density $R = 0.961$

- Total amount of fuel or liquid evaporated is proportional to heat release rate in combustion
- Important for combustion and spray drying

Test Matrix


Test Conditions	Ambient Pressure (kPa)	Fuel Temperature (°C)	Injection Pressure (MPa)
1	101	20	15
2	40	90	5
3	60	60	5
4	60	60	10
5	101	20	10
6	40	90	2

Injection at higher temperatures and lower than atmospheric pressures to simulate engine conditions. All tests carried out with calibrated gasoline-E10




Results

Six Plume Injector (20^o C)





1.035ms




Plume ID	Center R (mm)	Center θ (deg)	Center X (mm)	Center Y (mm)	Contour Area (mm ²)	Plume Angle (deg)	Peak Area (mm ⁻¹)	Total Area (mm ²)	% in Plume
1	12.1	19.8	8.0	5.1	33.2	7.4	0.69	11.3	13.2
2	11.8	77.5	-0.9	12.5	45.5	8.7	0.69	16.3	19.0
3	10.6	142.1	-11.8	7.5	43.6	8.5	0.69	14.7	17.1
4	10.3	218.0	-11.6	-5.3	47.4	8.9	0.75	15.4	17.9
5	12.6	278.0	-1.7	-11.5	44.8	8.6	0.71	16.4	19.0
6	11.8	330.8	6.8	-4.7	27.5	6.8	0.83	9.4	11.0
Total (plume separation at 15.9%)								83.5	97.2

Six Plume Injector (90° C)




- **Individual plumes not visible at low ambient pressures and high fuel temperatures**


High Speed Video



Six



Five





10 MPa, 20 °C, 101 KPa

2 MPa, 90 °C, 40 KPa



Five Plume Injector (20^o C)

1.035ms




Plume ID	Center R (mm)	Center θ (deg)	Center X (mm)	Center Y (mm)	Contour area (mm ²)	Plume angle (deg)	Peak area (mm ⁻¹)	Total area (mm ²)	% in plume
1	22.7	35.6	18.5	13.2	108.8	26.5	0.31	10.7	22.2
2	14.9	86.1	1.0	14.9	89.3	24.1	0.26	8.3	17.1
3	9.3	165.8	-9.0	2.3	103.4	25.8	0.32	9.3	19.3
4	13.6	253.4	-3.9	-13.0	99.8	25.4	0.28	8.6	17.8
5	21.3	305.6	12.4	-17.3	131.3	29.0	0.26	10.9	22.6
Total (plume separation at 5.3%)								47.8	99.1

Five Plume Injector (90° C)



- **Plumes still distinct at higher temperatures and lower ambient pressure**

Transient Total Surface Areas




- Profiles similar for same temperature and ambient pressures
- Lower injection pressures shows delay in peak values and lower peak values
- Injection temperature changes profile significantly

20 °C (1 and 5)

60 °C (3 and 4)


90 °C (2 and 6)

Comparison with Mechanical Patternator




Centriods from mechanical patternator is slightly further away

Comparison with Mechanical Patternator



Centroids from mechanical patternator is slightly further away



Conclusions

- **The peak surface area density increasing with increasing injection pressure.**
- **There is a merging of the plumes at higher temperatures and lower ambient pressures.**
- **The centroid locations obtained from the mechanical and optical patternator were in good agreement for both the injectors tested in this study**
- **The centroid locations do not vary much with pressure**
- **Extinction tomography is useful to characterizing the fuel spray from GDI injectors**